



Simplifying Alzheimer's Care

Ayelet Yehuda & Yaacov Levy-Dotan / Graphic Design

Infinite Boundaries



HomeFree™
at Home



HomeFree Systems Ltd.

Corporate Headquarters

2 Habarzel St.,
P.O.Box 13236, Tel Aviv 61132, Israel
Tel: +972-3-648 8690
Fax: +972-3-648 8691
E-mail: contact@homefreesys.com

US Office

P.O.Box 339
Yorktown Heights
NY 10598, USA
Tel: 1-800-606-0661
Fax: 914-962-1252

Visit us at: www.homefreesys.com



HomeFree™ at Home - Creating Safer Surroundings

For Alzheimer's patients and their caregivers, little is nearer to heart than safety, dignity and independence.

HomeFree at Home, an innovative wireless electronic system, is designed to alleviate concerns of wandering, by initiating alerts once a patient leaves a pre-defined area.

Simple, Silent Operation


HomeFree at Home is as simple as wearing a watch.

The wireless system comprises of:

- **Personal Watcher™ (PW)** - a uniquely designed watch, lightweight and waterproof
- **Home Monitoring Unit (HMU)** - a highly advanced data processing unit, featuring automatic communication and alarm capabilities
- **Caregivers' Pager** - for immediate local notification
- **Visual alerts** - optional strobe lights

The *PW* continuously sends a silent signal to the HMU, which processes the information and communicates with a central service at need.

When the patient leaves a designated area, the system automatically alerts caregivers. Alerts may be audible and/or visual and may also be sent locally by pager, where personal notification is preferred. Simultaneously, the HMU calls a central service, where a representative, in turn, contacts caregivers for additional assistance. Updates can also be sent via e-mail and Short Message Service (SMS) to cellular telephones, where available. Similar alerts are induced upon removal of the Personal Watcher.



Easiest Installation

HomeFree Personal Watcher is wrist worn, similarly to any other watch.

HomeFree at Home monitoring units are simply plugged into electricity and telephone outlets. No complicated installation required.



Central Monitoring Support

The **HomeFree at Home** solution is complemented by a sophisticated central monitoring system, enabling monitoring and support services. All local alarms are simultaneously received by a Central Monitoring Service.

Each patient is assigned a personalized schedule of daily routines and local range parameters, which can be automatically downloaded to the home unit or changed as needed.

The central service backs and supports home units by:

- updating caregivers' of patient's schedule.
- prompt verification by service representatives of alerts' reception by caregivers, assistance and guidance offered.
- confirmation of system's technical condition.



Full Compatibility with Care Centers Modules

HomeFree line of products includes a compatible monitoring system for care centers, **HomeFree at Care Centers**. Any *PW* wearer within a HomeFree monitored care center, is seamlessly transferred to the supervision of the center's system.

HomeFree Ltd. - a Trustworthy Business Partner

HomeFree is a leading provider of telehomecare systems. Company's systems are based on a unique technology, which design, development and production are ISO 9001 certified. **HomeFree** is a wholly owned subsidiary of the London Stock Exchange listed Dmatek Ltd. (DTK.L).